



Editorial Contact:
Patty Hoffman
Padilla Speer Beardsley
Tel: (612) 455-1781
Email: phoffman@padillaspeer.com



For a high resolution image visit:
www.colder.com/aseptiquikdc_photo

Colder Products Company Introduces AseptiQuik® DC Connectors

Single-Use Connection Technology for Both a Sterile Connection and Sterile Disconnection

St. Paul, Minn. – November 2, 2011 – [Colder Products Company](http://www.colder.com) announces the addition of the [AseptiQuik® DC Connectors](http://www.colder.com/aseptiquikdc) to its sterile connection offering for biopharmaceutical manufacturers seeking robust, sterile connection and disconnection for single-use processing. AseptiQuik DC Connectors – the first all-in-one single-use connection technology – deliver quick, easy and reliable connections for manufacturers using single-use systems and enable sterile media transfer, even in non-sterile environments.

The AseptiQuik DC's intuitive "CLICK-PULL-TWIST" design enables users to transfer media easily with less risk of operator error. After transfer is complete, the connector features a simple one-step disconnection that maintains media sterility by preventing external organisms from entering into the media flow path. The connector's robust design and automatic shutoff valves provide reliable performance without the need for clamps, fixtures or tube welders.

"Manufacturers have been requesting a sterile connection that can also be disconnected without risking contamination of the fluid path" said John Boehm, business unit manager, bioprocessing, Colder Products Company. "AseptiQuik DC connectors provide the security and reliability that users expect from single-use connection technology"

Colder's bioprocessing connections are manufactured in an ISO Class 7 certified cleanroom and meet USP Class VI and ADCF requirements. Contact Colder for validation test summaries, biocompatibility and extractables data.

Colder Products Company, a leader in single-use connection technology, will feature the AseptiQuik DC and other new products at BioProcess International (Booth #323) at the Long Beach Convention Center, California. Attendees have the opportunity to meet with our experts and discuss how to connect with Colder.

To learn more about how you can benefit from the AseptiQuik DC or other Colder bioprocessing connections please visit us online at www.colder.com/bio. For information on any of the other 7,500 connection solutions Colder offers, please visit us online at www.colder.com or call us at 800-444-2474.

###

(more)

About Colder Products Company

Colder Products Company is a leading provider of quick disconnect couplings and sterile connection technologies for the biopharmaceutical, medical and special industrial markets. For a free catalog CD or more information, contact: Colder Products Company, 1001 Westgate Drive, St. Paul, MN 55114. Phone: 651-645-0091. Toll-free: 800-444-2474. Fax: 651-645-5404. Web: www.colder.com.

CPC, Colder Products and Colder Products Company are registered trademarks with the U.S. Patent & Trademark office.